DCO Instructions for Use
Reusable Laparoscopic Ligation Clip Applier

EN – ENGLISH

IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference for surgical techniques.

DESCRIPTION

The GENICON Reusable Clip Appliers are designed to be used on single clips that have to be loaded individually. The overall length of the shaft of the instrument is approximately 12.99 inches (330mm). The clip applier is designed for introduction and use through an appropriately sized GENICON® trocar sleeve or larger sized sleeve.

INDICATIONS

The GENICON® reusable clip applier can be applied in endoscopic procedures for occlusion of vessels and other small tissue structures, as well as for radiographic markings.

CONTRAINDICATIONS

The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.

1. Do not use the GENICON® Clip Applier where metal ligating clips would not normally be used
2. This device is not intended for contraceptive tubal occlusion or renal artery ligation in nephrectomy. Device may be used to achieve hemostasis following transection of the fallopian tube.

WARNINGS AND PRECAUTIONS

1. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
2. The GENICON® reusable clip applier must be carefully inspected prior to and after each use. Devices with moving parts should be evaluated for proper function prior to application.
3. The reusable clip applier is provided nonsterile and must be cleaned and sterilized prior to each use.
4. Clip applier will arrive in open package. Do not use if package is damaged.
5. The clip applier and clips should not come in contact with electrosurgical electrodes.
6. The instrument is not intended for use as a grasper or dissector.
7. Appliers should only be used with GENICON clips.
8. Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to procedure.
9. Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
10. Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
11. Properly position the patient and note anatomical landmarks to introduced devices without unintended damage.
12. Do not use excessive force or in a manner not consistent with normal instrumentation use.
13. A thorough understanding of the principles and techniques involved in laparoscopic electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel, and damage to the device or other surgical instruments.
14. Verify that the devices are compatible with other products that will be used in surgery prior to the procedure.
15. All GENICON devices are latex free.
16. To avoid damage to clip cartridge, ensure the tips of applier are perpendicular to the cartridge when entering and withdrawing from the cartridge slot.
17. Detergents and solutions for reusable devices should have a PF between 7.0 and 9.5.
18. As a non-magnetic marking clip, it has low magnetic susceptibility and is safe for use in the MRI environment. Positional errors and artifacts may increase dramatically with increasing static magnetic field strength and/or increasing field of view size. Artifact size may vary greatly for different MRI pulse sequences.
19. The titanium clip is biologically inert and can remain in body without consequences.

To prevent premature closure or release of the clip:

- Never grasp the applier by both shanks when loading the clip.
- Do not squeeze applier handles before jaws are placed around vessel or tissue.
- Ensure appliers are fully closed to handle stop point when applying clip.
- Inspect ligation site to ensure secure application and hemostasis.

To prevent damage to applier jaws:

- Avoid introducing applier tips at an angle into the cartridge rack, this may misalign jaws.
- Inspect applier and jaws before each use.
- Use only GENICON branded clips.
- Open applier for thorough cleaning and sterilization.

CONTENTS

P/N	DESCRIPTION	SIZE	CLIP SIZE
350-005-000	Single-Action Jaw Non-Detachable Ligation Clip Applier	5 x 330 mm	Small
350-010-000	Double-Action Jaw Detachable Ligation Clip Applier	10 x 330 mm	Large
350-010-002	Double-Action Jaw Detachable Ligation Clip Applier	10 x 330 mm	Medium/Large
350-010-003	Single-Action Jaw Detachable Ligation Clip Applier	10 x 330 mm	Medium/Large
350-010-090	Single-Action Jaw Non-Detachable Ligation Clip Applier	10 x 330 mm	Medium

INSTRUCTIONS


1. Make sure the device has been sterilized before each use including the first use.
2. Select applier according to clip size.

3. Place clip cartridge firmly on a stable surface and hold clip cartridge in place while removing clips.
4. Grasp the applier by the handle.
5. Make sure the tips of applier are perpendicular to the cartridge.
6. Insert the applier over the clip until it stops. Do not force the applier. It should enter and withdraw from the cartridge slot smoothly.
7. Retract the applier away from the cartridge. The clip should be secure in the applier jaws. It is not necessary to maintain trigger tension to hold the clip in the applier jaws.
8. Insert the clip applier into the suitable trocar sleeve.
9. Place the jaws around the tissue to be ligated. Rotate the applier shaft to facilitate placement of the clip. Ensure that the tissue to be ligated is located completely within the confines of the clip prior to closure.
10. Squeeze the handle firmly until it will not move further without excessive force. Failure to fully squeeze the handle may result in malformed clips which may result in incomplete closure and lack of hemostasis.
11. Release the handle and remove applier. If handle is not fully decompressed clip may remain on applier resulting in improper ligation, injury or other serious surgical complications.
12. After firing the clip, retract the applier from the trocar sleeve.

CLEANING INSTRUCTIONS

Detachable clip applier must be disassembled (two pieces – applier and insert) and thoroughly cleaned and sterilized prior to EACH use. This includes first use. Non-Detachable clip appliers are not disassembled for sterilization.

DISASSEMBLY

The clip applier disassembles into two pieces. Firmly grasp the handle in your left hand while placing the thumb and forefinger of your right hand around the prongs. Turn your right hand counter-clockwise (towards body) to release insert **1**. After loosening the insert, pull slightly to remove insert from applier **2**. If when you pull, the insert does not move; continue turning the prongs counter-clockwise. When complete, there will be two parts; the applier and the insert **3**.

CLEANING

The user must ensure that cleaning and sterilization are conducted in accordance with the appropriate guidelines, standards or national Health Authorities' requirements.

Applier jaws are delicate and can easily become damaged. Mishandling of appliers may result in improper load and/or closure of the ligating clips. Appropriate care, cleaning and maintenance are important to ensure proper function. All GENICON's reusable clip appliers are supplied non-sterile. The instrument should be cleaned and sterilized prior to each use. Clean the instrument before sterilization and in the same manner as any reusable instrument and in accordance with hospital practice.

- We recommend using the ultrasonic cleaning method to clean the reusable clip applier.
- It is recommended to be disinfected by 2% alkaline glutaraldehyde. Please follow the instructions of the disinfectant manufacturer.
- Dry the shaft channel by compressed air before steam sterilization. We recommend using the steam sterilization.


Sterilization Method	Product Category	Temperature	Minimum Time	Pressure
Gravity Cycle	Appliance	121°C	30 Min	102.9 Kpa
Pre Vacuum Steam Cycle	Appliance	132~134°C	4 Min	205.8 Kpa

NOTE: The Sterilizer Manufacturer's Instructions for operations and load configuration should be followed explicitly.

NOTE: The clip applier may be autoclaved up to 200 times using the above procedure.

ASSEMBLY

This applies to the detachable clip applicator. The non-detachable clip applicator does not require assembly.

The detachable clip applicator is composed of two pieces, the applicator (A) and the insert (B). Firmly grasp the applicator handle in your left hand and the insert in your right hand **1**. Line the insert up with the opening of the applicator. Gently push the insert into the applicator until you feel resistance **2**. At that point, begin twisting the insert in a clockwise manner – away from your body **3**. Continue twisting the insert until it is snug (it is locked when it can no longer twist).

STORAGE

- Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65 % (non-condensing). Handle with care.
- Keep safe in clean environment at max 80% RH with good ventilation, which is kept away from the sun, fire, and toxic air.
- The clip applicator should be stored assembled to avoid misplacing the insert.

SHELF LIFE

- Do not use if package is damaged.
- The manufacturer declines responsibility for problems resulting from improper use.
- Prices, specifications, and model availability are subject to change without notice.

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending. GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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© GENICON® Applier_AD Page 5 of 30