Instructions for Use

Titanium Hemostatic Clip

EN – ENGLISH

BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or re-sterilization of this device may create risk of contamination and patient infection.
5. Do not reuse, reprocess or re-sterilize this device

DESCRIPTION

The titanium hemostatic clip is supplied packaged EO sterile and has been specifically designed to securely grip the vessel or tissue to provide complete atraumatic occlusion. User shall select appropriate clip size according to procedure ensuring that all tissues are completely encompassed by clip.

INDICATIONS

The GENICON implantable clips are intended for use during procedures for the purpose of approximating soft tissue, closing off vessels and other structures in order to stop bleeding or to connect internal tissues to aid in healing.

CONTRAINDICATIONS

The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.

WARNINGS AND PRECAUTIONS

1. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
2. GENICON single-use products are intended for single patient use only – DO NOT RESTERILIZE.
3. Do not use any pre-sterilized products if the packaging sterile barrier is damaged
4. The clip applicator and clips should not come in contact with electrosurgical electrodes.
5. Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to procedure.
6. Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
7. Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
8. Properly position the patient and note anatomical landmarks to introduced devices without unintended damage.
9. Do not use excessive force or in a manner not consistent with normal instrumentation use.
10. A thorough understanding of the principles and techniques involved in laparoscopic electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel, and damage to the device or other surgical instruments.
11. Verify that the devices are compatible with other products that will be used in surgery prior to the procedure.
12. All GENICON devices are latex free.
13. Dispose of all used or damaged products using standard hospital practices for biohazard control.
14. To avoid damage to clip cartridge, ensure the tips of applicator are perpendicular to the cartridge when entering and withdrawing from the cartridge slot.
15. As a non-magnetic marking clip, it has low magnetic susceptibility and is safe for use in the MRI environment. Positional errors, phantoms, and artifacts may increase dramatically with increasing static magnetic field strength and/or increasing field of view size. Artifact size may vary greatly for different MRI pulse sequences.
16. The titanium clip is biologically inert and can remain in body without consequences.
17. Do not use the GENICON® Clip Applicator where metal ligating clips would not normally be used
18. This device is not intended for contraceptive tubal occlusion or renal artery ligation in nephrectomy. Device may be used to achieve hemostasis following transection of the fallopian tube.
19. Federal (USA) law restricts this device to sale by or on the order of a physician.

To prevent premature closure or release of the clip:

- Never grasp the applicator by both shanks when loading the clip.
- Do not squeeze applicator handles before jaws are placed around vessel or tissue.
- Ensure applicators are fully closed to handle stop point when applying clip.
- Inspect ligation site to ensure secure application and hemostasis.
- Do not apply a staple or staple line over a clip

To prevent damage to applicator jaws:

- Avoid introducing applicator tips at an angle into the cartridge rack, this may misalign jaws.
- Inspect applicator and jaws before each use.
- Use only GENICON branded clips.
- Use recommended applicators.
- Open applicator for thorough cleaning and sterilization.
- This product should not be used on patients with Nickel or Titanium allergies

CONTENT

| P/N | DESCRIPTION | COLOR CODE |
|---|---------------------------------------|------------|
| 350-000-101 use with 350-005-000, 330-005-000 | 20 Cartridges of 6 Small Clips | Blue |
| 350-000-102 use with 350-010-090, 330-005-002, 330-005-003 | 20 Cartridges of 6 Medium Clips | White |
| 350-000-103 use with 350-010-002, 350-010-003, 330-005-004 | 20 Cartridges of 6 Medium/Large Clips | Green |
| 350-000-104 use with 350-010-000, 350-010-001, 330-005-005 | 20 Cartridges of 6 Extra Large Clips | Yellow |

INSTRUCTIONS

1. Select applier according to clip size.
2. Place clip cartridge firmly on a stable surface and hold clip cartridge in place while removing clips.
3. Grasp the applier by the handle.
4. Make sure the tips of applier are perpendicular to the cartridge.
5. Insert the applier over the clip until it stops. Do not force the applier. It should enter and withdraw from the cartridge slot smoothly.
6. Retract the applier away from the cartridge. The clip should be secure in the applier jaws. It is not necessary to maintain trigger tension to hold the clip in the applier jaws.
7. Insert the clip applier into the suitable trocar sleeve.
8. Place the jaws around the tissue to be ligated. Rotate the applier shaft to facilitate placement of the clip. Ensure that the tissue to be ligated is located completely within the confines of the clip prior to closure.
9. Squeeze the handle firmly until it will not move further without excessive force. Failure to fully squeeze the handle may result in malformed clips which may result in incomplete closure and lack of hemostasis.
10. Release the handle and remove applier. If handle is not fully decompressed clip may remain on applier resulting in improper ligation, injury or other serious surgical complications.
11. After firing the clip, retract the applier from the trocar sleeve.

STORAGE

Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65 % (non-condensing). Handle with care.

Before using the device, please check the parts to be inserted into human body to make sure the surface is smooth, without sharp edges or burs which possibly cause safety damage.

SHELF LIFE

- Do not use if package is damaged.
- The ligating clips are intended for single use and cannot be resterilized.
- The manufacturer declines responsibility for problems resulting from improper use.
- Prices, specifications, and model availability are subject to change without notice.

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.
GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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