




## GENICON SINGLE USE MONOPOLAR INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Monopolar Ratcheting and Non-Ratcheting Dissectors, Scissors and Graspers

**BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

**IMPORTANT!**

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or resterilization of this device may create risk of contamination and patient infection.
5. Do not reuse, reprocess or resterilize this device.

**DESCRIPTION**

The GENICON Monopolar Laparoscopic Instruments are sterile packaged Single Use Mono-Polar instruments including graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue.

**INDICATIONS FOR USE –**

Endoscopic surgical procedures. It is a family of instruments which includes graspers, dissectors, and scissors, which are intended to be used to grasp, manipulate, cut and cauterize soft tissue.

**CONTRAINDICATIONS**

1. The GENICON monopolar Laparoscopic single use instruments are NOT intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transaction of the fallopian tube.
2. These devices are intended for use only as indicated.

**WARNINGS AND PRECAUTIONS**

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.
2. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
3. A thorough understanding of the principles and techniques involved in electrosurgical procedures is necessary to avoid shock and burn to both patient and operator. Verify compatibility of instrumentation, and ensure that electrical insulation or grounding is not compromised.
4. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
5. Do not use any pre-sterilized products if the packaging sterile barrier is damaged.
6. Utilize available technology (tissue response generator, active electrode monitoring) to further eliminate concerns about insulation failure, capacitive coupling and interference with other electrical equipment.
7. Cords should not be bundled, wrapped around metal instruments affixed to the patient or placed directly on patient's skin. Antenna coupling may occur when running the electrosurgical cord parallel to other patient or user contact cords or leads.
8. Do not place electrocautery cable on camera cable to avoid monitor display interference.
9. Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use. Do not use if inspection or current leakage tests indicate damage. Insulation failures may result in burns or other injuries to the patient or operator.
10. Single Use devices may not be repaired, modified or re-processed.
11. Introduce the instrument through the cannula carefully to avoid damaging the working tip.
12. Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65 % (non-condensing). Handle with care.
13. DO NOT USE in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N<sub>2</sub>O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
14. Adjust the HF-output power to match with the intended procedure. Use the lowest possible power setting that will achieve the desired surgical effect. Such as the low-voltage waveform (cut) and brief intermittent activation. This will reduce the potential for capacitive coupling and/or inadvertent burning of tissues.
15. Do not use power settings that may result in more than 120W being delivered to the GENICON monopolar laparoscopic instruments. Excessive power levels may result in instrument malfunction and possible patient or user injury. The rated voltage for these device is 5000 Vp.
16. Compatibility of parts from different manufacturers must be verified prior to conducting the procedure.
17. The devices are intended for use with electrosurgical generators and accessories in compliance with safety standards IEC 60601-1, IEC 60601-1-2, and IEC 60601-2-2 including applicable national/regional differences. Refer to the electrosurgical generator manual to verify compatibility and corresponding settings. Follow all safety precautions.
18. Ensure device grasping or cutting surfaces are fully visible prior to engaging the electrical current to avoid unintended results. Keep the working-end under full and unobstructed visualization during use.
19. Damage to instrument may occur if attempting to cut staples, clips, or other non tissue based materials.
20. Use appropriate technique to achieve hemostasis if not present after removal of instrument.
21. Do not activate in close proximity or direct contact with another instrument, such as clips, staple, laparoscope, etc. Do not active the instruments in contact with, in close proximity to, or in conjunction with conductive fluids (e.g. blood or saline). This may alter the electrical path and cause unintended tissue damage.
22. Electrosurgical devices used in conjunction with a laser/Argon beam may create the potential for the development of a gas embolism.
23. Do not activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling with other surgical equipment or cause damage to the device
24. Electrosurgical generators may cause unintended destruction of tissue and are dangerous if operated improperly. Follow all instructions for use required by the generator manufacturer.
25. Mono-polar products should be connected ONLY to a mono-polar power connection on the generator.
26. Do not use excessive force or in a manner not consistent with normal instrumentation use.
27. Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.
28. The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.
29. Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
30. Keep the active electrodes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning. Injury to operating room personnel may result.
31. When the active electrode is not being used, it should be placed in an insulated holster. Do not place instrument on the patient or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire or patient burns.

**POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS**

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

**INSTRUCTIONS FOR USE**

1. Open the package by using standard sterile technique.
2. Inspect the unit for defects – do not use if unit is defective. **Refer to picture 1**
  - A. Shaft
  - B. Color-Coded Knob
  - C. Monopolar Post
  - D. Ratchet (*ratcheting version only*)
  - E. Color-Coded Thumb Loop
  - F. Handle
3. When attaching a mono-polar generator cable to the cautery connector located on the top of the instrument **use an electro-surgical cable that has a 4mm female connector** and follow the generator manufacturer's instructions for use and setup. **Refer to picture 2**
4. To rotate the shaft, turn the knob either direction to desired location. **Refer to picture 3**
5. For ratcheting instruments, position tissue inside the jaws and squeeze handle to desired clamping tightness. To activate the ratchet rotate the ratchet arm into the down position. **Refer to picture 4**. To deactivate the ratchet rotate the ratchet arm into the up position Spread handle to release the tissue. **Refer to picture 5**
6. Keep the contact surface of the instrument clean during the operation. Wipe off any dried residue.
7. The handle and jaw mechanisms are designed to open smoothly. Care should be taken not to forcibly open the jaws wider as undue stress may damage the handle, insulation, or jaw components.

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Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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