







## GENICON SINGLE USE MONOPOLAR INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Monopolar Ratcheting and Non-Ratcheting Dissectors, Scissors and Graspers

## **BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

### **IMPORTANT!**

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or reesterilization of this device may create risk of contamination and patient infection.
5. Do not reuse, reprocess or reesterilize this device.

### **DESCRIPTION**

The GENICON Monopolar Laparoscopic Instruments are sterile packaged Single Use Mono-Polar instruments including graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue.

### **INDICATIONS FOR USE –**

Endoscopic surgical procedures. It is a family of instruments which includes graspers, dissectors, and scissors, which are intended to be used to grasp, manipulate, cut and cauterize soft tissue.

### **CONTRAINDICATIONS**

1. The GENICON monopolar Laparoscopic single use instruments are NOT intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transaction of the fallopian tube.
2. These devices are intended for use only as indicated.

### **WARNINGS AND PRECAUTIONS**

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.
2. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
3. A thorough understanding of the principles and techniques involved in electrosurgical procedures is necessary to avoid shock and burn to both patient and operator. Verify compatibility of instrumentation, and ensure that electrical insulation or grounding is not compromised.
4. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
5. Do not use any pre-sterilized products if the packaging sterile barrier is damaged.
6. Select a generator with necessary safety mechanisms to eliminate interference with other electrical equipment.
7. Do not loop electrocautery cable or place directly on patient's skin.
8. Do not place electrocautery cable on camera cable to avoid monitor display interference.
9. Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use. Do not use if inspection or current leakage tests indicate damage.
10. Single Use devices may not be repaired, modified or re-processed.
11. Introduce the instrument through the cannula carefully to avoid damaging the working tip.
12. Store all products safely in a climate controlled environment and handle with care.
13. A risk of injury may arise from any sparks or ignition of combustible gases. Be sure to follow the safety information in the operating instructions of the HF generator.
14. Adjust the HF-output power to match with the intended procedure. Voltage/power levels should be set as low as possible to achieve the desired effects. This will reduce the potential for capacitive coupling and/or inadvertent burning of tissues.
15. Do not use power settings that may result in more than 2600 Volts peak being delivered to the GENICON monopolar laparoscopic instruments. Excessive power levels may result in instrument malfunction and possible patient or user injury.
16. Compatibility of parts from different manufacturers must be verified prior to conducting the procedure.
17. Ensure device grasping or cutting surfaces are fully visible prior to engaging the electrical current to avoid unintended results. Keep the working-end under full and unobstructed visualization during use.
18. Damage to instrument may occur if attempting to cut staples, clips, or other non tissue based materials.
19. Use appropriate technique to achieve hemostasis if not present after removal of instrument.
20. Ensure that the energized tip is not in contact with a conductive irrigation fluid, staples, clips, or other conductive device. Users are cautioned that activating the mono-polar device simultaneously with suction/irrigation devices may alter the path of the electrical energy away from the target tissue.
21. Mono-polar devices used in conjunction with a laser/Argon beam may create the potential for the development of a gas embolism.
22. Capacitive coupling may occur if the device is activated and not in position to deliver energy to the target tissue.
23. Electrosurgical generators may cause unintended destruction of tissue and are dangerous if operated improperly. Follow all instructions for use required by the generator manufacturer.
24. Mono-polar products should be connected ONLY to a mono-polar power connection on the generator.

## POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

## INSTRUCTIONS FOR USE

### Common/General Instructions for Use

1. Open the package by using standard sterile technique.
2. Remove protective tip from the distal end of the assembly and the monopolar post protector from the handle. **Refer to pictures 1 and 2**
3. Inspect the unit for defects – do not use if unit is defective
  - 3.1. **Refer to picture 3**
    - A. Shaft
    - B. Color-Coded Knob
    - C. Monopolar Post
    - D. Handle
  - 3.2. **Refer to picture 4**
    - A. Shaft
    - B. Color-Coded Knob
    - C. Monopolar Post
    - D. Handle
    - E. Ratchet
4. Ensure that the product functions and is intact. When attaching a mono-polar generator cable to the cautery connector located on the top of the instrument follow the generator manufacturer's instructions for use and setup. **Refer to picture 5**
5. To rotate the shaft, turn the knob either direction to desired location. **Refer to picture 6**
6. For ratcheting instruments, position tissue inside the jaws and squeeze handle to desired clamping tightness. Press the release button in front of the handle to release. Spread handle to release the tissue. **Refer to pictures 7 and 8**
7. Keep the contact surface of the instrument clean during the operation. Wipe off any dried residue.
8. The handle and jaw mechanisms are designed to open smoothly. Care should be taken not to forcibly open the jaws wider as undue stress may damage the handle, insulation, or jaw components.

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.

GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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