













## SINGLE USE SIMPLE RETRIEVAL BAG INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Simple Retrieval bag

## **BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

### **IMPORTANT!**

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or re-sterilization of this device may create risk of contamination and patient infection.
5. Do not reuse, reprocess or re-sterilize this device

### **DESCRIPTION**

The GENICON Single-Use Simple Retrieval Bag is comprised of a flexible plastic bag with a large, easily accessible opening, a top and bottom tab, and a closure suture. A string with a closure suture facilitates closure of the specimen bag after the specimen had been collected. This device is packaged and sterilized for single use only. Do not re-use, reprocess, or re-sterilize. Discard after use.

### **INDICATIONS FOR USE**

The GENICON Single-Use Simple Retrieval Bag is indicated for use in endoscopic procedures to allow tissues or debris to be removed from the abdominal cavity. These devices have applications in gynecologic, general, thoracic, and urologic endoscopic procedures.

### **CONTRAINDICATIONS**

- The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.

### **WARNINGS AND PRECAUTIONS**

1. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
2. GENICON single-use products are intended for single patient use only – DO NOT RESTERILIZE.
3. Do not use any pre-sterilized products if the packaging sterile barrier is damaged
4. Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to procedure.
5. Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
6. Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
7. Properly position the patient and note anatomical landmarks to introduced devices without unintended damage.
8. Do not use morcellators with this device
9. Care should be taken to avoid contact of the bag with sharp instruments, cutting devices, and electrosurgical and laser instrument.
10. This device is not intended for use with any tissue that will not fit within the confines of the specimen bag and allow complete closure of the bag.
11. Do not attempt to remove the specimen through the trocar or cannula as this may lead to bag rupture and spillage of contents.
12. The specimen bag includes a radiopaque thread that is clearly visible on CT Scans and may be visible in standard, plain radiography.
13. Excessive forces should be avoided during bag extraction. If the bag and its contents are too large to be extracted, carefully enlarge the access site to facilitate easy bag removal.
14. If the bag and its contents are too large to be extracted, carefully enlarge the access site to facilitate easy bag removal.
15. Verify that the devices are compatible with other products that will be used in surgery prior to the procedure.
16. Store device safely in a climate controlled environment at a temperature less than 54° C and relative humidity less than 80%.
17. The GENICON Simple Retrieval Bag is latex free and sterilized using EtO sterilization, except as noted.
18. Dispose of all used or damaged products using standard hospital practices for biohazard control.
19. Do not introduce The GENICON Simple Retrieval Bag through cannula smaller than 10 mm.

### **POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS**

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

### **INSTRUCTIONS FOR USE – 550-000-009 | 550-000-010 | 550-000-017**

1. Inspect all components carefully for damage during shipment, handling or after maintenance. **Refer to picture 1a**
  - A. Closure Suture
  - B. Top tab
  - C. Radiopaque thread
  - D. Specimen bag
2. Inspect sterile packaging, and verify that sterile barriers are intact. Do not use if package is compromised.
3. Verify that the sizes of all surgical components selected are compatible.

4. Open package using sterile technique and place instrument on sterile field.
5. Grasp the top tab with an atraumatic laparoscopic grasper and insert the simple bag completely through the desired 10 mm or larger port; applying sterile saline solution to the bag will aid in introduction **Refer to picture 2a**
6. Place specimen into the simple bag. **Refer to picture 3a**
7. Once the specimen is completely within the bag, pull the closure suture with a second atraumatic laparoscopic grasper until the bag is completely cinched closed. **Refer to picture 4a**
8. Once extraction is desired, prepare by pulling the specimen bag with an atraumatic laparoscopic grasper until the leading edge of the bag is inside the tip of the cannula. **Refer to picture 5a**
9. While maintaining the grip on the bag, pull the cannula and instrument out of the port site together until the top of the closed bag reaches the body surface or access site. Carefully remove the bag from the access site. **Refer to picture 6a**

**NOTE:** Do not attempt to pull the specimen through the cannula.  
**NOTE:** If the bag and its contents are too large to be extracted, carefully enlarge the access site to facilitate easy bag removal.

#### INSTRUCTIONS FOR USE – 550-000-011 | 550-000-012 | 550-000-013

1. Inspect all components carefully for damage during shipment, handling or after maintenance. **Refer to picture 1b**
  - A. Closure Suture
  - B. Top tab
  - C. Radiopaque thread
  - D. Specimen bag
  - E. Bottom tab
2. Inspect sterile packaging, and verify that sterile barriers are intact. Do not use if package is compromised.
3. Verify that the sizes of all surgical components selected are compatible.
4. Open package using sterile technique and place instrument on sterile field.
5. Grasp the bottom tab (E) with an atraumatic laparoscopic grasper and roll the short side of the bag along the shaft of the grasper. **Refer to picture 2b**
6. Insert the rolled bag into the appropriate size cannula as below **Refer to picture 3b**

Product Code	Compatible GENICON Cannula
550-000-011	12 mm or larger
550-000-012	15 mm
550-000-013	15 mm
7. Once the bag material is visible inside the cavity, using a second grasper grab the available tab (Top Tab A) and release the first grasper to pull the bag into the cavity. **Refer to picture 4b**
8. Place specimen into the simple bag. **Refer to picture 5b**
9. Once the specimen is completely within the bag, pull the closure suture with a second atraumatic laparoscopic grasper until the bag is completely cinched closed. **Refer to picture 6b**
10. Once extraction is desired, prepare by pulling the specimen bag with an atraumatic laparoscopic grasper until the leading edge of the bag is inside the tip of the cannula. **Refer to picture 7b**
11. While maintaining the grip on the bag, pull the cannula and instrument out of the port site together until the top of the closed bag reaches the body surface or access site. Carefully remove the bag from the access site. **Refer to picture 8b**

**NOTE:** Do not attempt to pull the specimen through the cannula.  
**NOTE:** If the bag and its contents are too large to be extracted, carefully enlarge the access site to facilitate easy bag removal.

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.  
GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

**Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products:** Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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