
## GENICON SINGLE USE MONOPOLAR X-SURGE SHAFT INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Monopolar Dissector, Scissor and Grasper Shafts. Handle not included.

**BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

### IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques
2. This device was designed, tested and manufactured for single use patient use only
3. Reuse of this device may lead to its failure and subsequent patient injury
4. Reprocessing and/or resterilization of this device may create risk of contamination and patient infection
5. Do not reuse, reprocess or resterilize this device

### DESCRIPTION

The GENICON X-Surge Monopolar Laparoscopic Instruments are sterile packaged single-use monopolar attachments intended for use in combination with a compatible handle (*supplied separately*). It includes graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue (*handle required*).

### INDICATIONS FOR USE-

Endoscopic surgical procedures, it is a family of instruments which include graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue.

### CONTRAINDICATIONS

1. The GENICON monopolar Laparoscopic X-Surge single use instruments are NOT intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transaction of the fallopian tube.
2. These devices are intended for use only as indicated.

### WARNINGS AND PRECAUTIONS

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for The use of GENICON® laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON® products are noted in the specific sections.
2. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure. Consult medical literature or country specific regulations for specific techniques, complications, and hazards prior to procedure.
3. A thorough understanding of the principles and techniques involved in electrosurgical procedures is necessary to avoid shock and burn to both patient and operator. Verify compatibility of instrumentation, and ensure that electrical insulation or grounding is not compromised.
4. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
5. Do not use any pre-sterilized products if the packaging sterile barrier is damaged. Handle with care.
6. Select a generator with necessary safety mechanisms to eliminate interference with other electrical equipment.
7. Do not loop electrocautery cable or place directly on patient's skin.
8. Do not place electrocautery cable on camera cable to avoid monitor display interference.
9. The device must be carefully inspected prior to and after each use, including compatible handles and related accessories. Devices with moving parts should be evaluated for proper function prior to application. Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use. Do not use if inspection or current leakage tests indicate damage.
10. Single Use devices may not be repaired, modified or re-processed.
11. Properly position the patient and note anatomical landmarks to introduced devices without unintended damage. Introduce the instrument through the cannula carefully to avoid damaging the working tip.
12. Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65 % (non-condensing). Handle with care.
13. A risk of injury may arise from any sparks or ignition of combustible gases. Be sure to follow the safety information in the operating instructions of the HF generator.
14. Adjust the HF-output power to match with the intended procedure. Voltage/power levels should be set as low as possible to achieve the desired effects. This will reduce the potential for capacitive coupling and/or inadvertent burning of tissues.
15. Do not use power settings that may result in more than 120W being delivered to the devices. Excessive power levels may result in instrument malfunction and possible patient or user injury. The rated voltage for these devices is 3800 Vp.
16. Compatibility of parts from different manufacturers must be verified prior to conducting the procedure.
17. The devices are intended for use with electrosurgical generators and accessories in compliance with safety standards IEC 60601-1, IEC 60601-1-2, and IEC 60601-2-2 including applicable national/regional differences. Refer to the electrosurgical generator manual to verify compatibility and corresponding settings. Follow all safety precautions.
18. Ensure device cutting and coagulation surfaces are fully visible prior to engaging the electrical current to avoid unintended results. Keep the working-end under full and unobstructed visualization during use.
19. Use appropriate technique to achieve hemostasis if not present after removal of instrument.
20. Ensure that the energized tip is not in contact with a conductive irrigation fluid, staples, clips, or other conductive device. Users are cautioned that activating the monopolar function simultaneously with suction/irrigation function may alter the path of the electrical energy away from the target tissue.
21. Monopolar devices used in conjunction with a laser/Argon beam may create the potential for the development of a gas embolism.
22. Capacitive coupling may occur if the device is activated and not in position to deliver energy to the target tissue.
23. Electrosurgical generators may cause unintended destruction of tissue and are dangerous if operated improperly. Follow all instructions for use required by the generator manufacturer.
24. Monopolar products should be connected ONLY to a monopolar power connection on the generator.
25. Care must be taken when using laparoscopic instrumentation to avoid damage to major vessels and other anatomic structures.
26. Establish and maintain adequate pneumoperitoneum to reduce the risk of injury to internal structures.
27. Do not use excessive force or in a manner not consistent with normal instrumentation use.
28. Cardiac Pacemakers may be damaged by electrosurgical current. Never use electrosurgical current on patients fitted with pacemakers. Consult a cardiologist for intervention.

### POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

### INSTRUCTIONS FOR USE

1. Inspect all components carefully for any damage that may have occurred during shipment. Inspect packaging marked as sterile, and verify that sterile barriers are intact. Do not use if sterile package is damaged or opened.
2. Open the package using standard sterile technique

3. Inspect the unit for defects – do not use if unit is defective. **Refer to picture 1**
  - A. Jaws
  - B. Shaft/Insulation
  - C. Locking Knob
  - D. Drive Rod/Ball Joint
4. Open the handle to fully extend the attachment socket. **Refer to picture 2**
5. Align the Ball Joint of the shaft and insert into the exposed attachment socket of the handle. The shaft is to be parallel to the working length of the handle. **Refer to picture 3**
6. Close handle to temporarily secure the shaft in place, and slide the Locking Knob to the proximal end of the shaft. **Refer to picture 4**
7. Rotate the Locking Knob clockwise until tightened to handle, ensure fit is secure. **Refer to picture 5**
8. Using the handle, open the tip to full aperture and verify proper assembly and working condition. **Refer to picture 6**
9. Keep the contact surface of the instrument clean during the operation. Wipe off any dried residue.
10. The tip mechanism is designed to open and close smoothly. Care should be taken not to forcibly open the jaws wider as undue stress may damage the insulation, or jaw components.

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Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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