




GENICON SINGLE-USE MONOPOLAR X-SURGE TIPS INSTRUCTIONS FOR USE

Including sterile packaged Single-Use Monopolar Dissector, Scissor and Grasper Tips. Handle not included.

BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques
2. This device was designed, tested and manufactured for single use patient use only
3. Reuse of this device may lead to its failure and subsequent patient injury
4. Reprocessing and/or reesterilization of this device may create risk of contamination and patient infection
5. Do not reuse, reprocess or reesterilize this device

DESCRIPTION

The GENICON X-Surge Monopolar Laparoscopic Instruments are sterile packaged single-use monopolar attachments intended for use in combination with a compatible handle (*supplied separately*). It includes graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue (*handle required*).

INDICATIONS FOR USE

Endoscopic surgical procedures, it is a family of instruments which include graspers, dissectors, and scissors, intended to grasp, manipulate, cut, and cauterize soft tissue.

CONTRAINDICATIONS

1. The GENICON monopolar Laparoscopic X-Surge single use instruments are NOT intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transaction of the fallopian tube.
2. These devices are intended for use only as indicated.

WARNINGS AND PRECAUTIONS

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.
2. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
3. A thorough understanding of the principles and techniques involved in electrosurgical procedures is necessary to avoid shock and burn to both patient and operator. Verify compatibility of instrumentation, and ensure that electrical insulation or grounding is not compromised.
4. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
5. Do not use any pre-sterilized products if the packaging sterile barrier is damaged.
6. Utilize available technology (tissue response generator, active electrode monitoring) to further eliminate concerns about insulation failure, capacitive coupling and interference with other electrical equipment.
7. Cords should not be bundled, wrapped around metal instruments affixed to the patient or placed directly on patient's skin. Antenna coupling may occur when running the electrosurgical cord parallel to other patient or user contact cords or leads.
8. Do not place electrocautery cable on camera cable to avoid monitor display interference.
9. Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use. Do not use if inspection or current leakage tests indicate damage. Insulation failures may result in burns or other injuries to the patient or operator.
10. Single Use devices may not be repaired, modified or re-processed.
11. Introduce the instrument through the cannula carefully to avoid damaging the working tip.
12. Transport and store all products safely at ambient temperatures of 5°C to 30°C and relative humidity of 35% to 65 % (non-condensing). Handle with care.
13. DO NOT USE in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
14. Adjust the HF-output power to match with the intended procedure. Use the lowest possible power setting that will achieve the desired surgical effect. Such as the low-voltage waveform (cut) and brief intermittent activation. This will reduce the potential for capacitive coupling and/or inadvertent burning of tissues.
15. Do not use power settings that may result in more than 120W being delivered to the GENICON monopolar laparoscopic instruments. Excessive power levels may result in instrument malfunction and possible patient or user injury. The rated voltage for these device is 3800 Vp.
16. Compatibility of parts from different manufacturers must be verified prior to conducting the procedure.
17. The devices are intended for use with electrosurgical generators and accessories in compliance with safety standards IEC 60601-1, IEC 60601-1-2, and IEC 60601-2-2 including applicable national/regional differences. Refer to the electrosurgical generator manual to verify compatibility and corresponding settings. Follow all safety precautions.
18. Ensure device grasping or cutting surfaces are fully visible prior to engaging the electrical current to avoid unintended results. Keep the working-end under full and unobstructed visualization during use.
19. Damage to instrument may occur if attempting to cut staples, clips, or other non tissue based materials.
20. Use appropriate technique to achieve hemostasis if not present after removal of instrument.
21. Do not activate in close proximity or direct contact with another instrument, such as clips, staple, laparoscope, etc. Do not active the instruments in contact with, in close proximity to, or in conjunction with conductive fluids (e.g. blood or saline). This may alter the electrical path and cause unintended tissue damage.
22. Electrosurgical devices used in conjunction with a laser/Argon beam may create the potential for the development of a gas embolism.
23. Do not activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling with other surgical equipment or cause damage to the device
24. Electrosurgical generators may cause unintended destruction of tissue and are dangerous if operated improperly. Follow all instructions for use required by the generator manufacturer.
25. Mono-polar products should be connected ONLY to a mono-polar power connection on the generator.
26. Do not use excessive force or in a manner not consistent with normal instrumentation use.
27. Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.
28. The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.
29. Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
30. Keep the active electrodes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning. Injury to operating room personnel may result.
31. When the active electrode is not being used, it should be placed in an insulated holster. Do not place instrument on the patient or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire or patient burns.
32. The device must be carefully inspected prior to and after each use, including compatible handles and related accessories. Devices with moving parts should be evaluated for proper function prior to application. Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use. Specific attention is called to the O-ring on the distal end of the handle assembly. Do not use if inspection or current leakage tests indicate damage.
33. Failure to completely close the jaws and/or handle prior to assembly may compromise the cutting/grasping performance.

34. No gap should remain between the tip and shaft after assembly. Failure to properly assemble the device may result in instrument malfunction and possible patient or user injury.

POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

INSTRUCTIONS FOR USE

1. Inspect all components carefully for any damage that may have occurred during shipment. Inspect packaging marked as sterile, and verify that sterile barriers are intact. Do not use if sterile package is damaged or opened.
2. Open the package using standard sterile technique
3. Inspect the unit for defects – do not use if unit is defective. **Refer to picture 1**
 - A. Tip Protector
 - B. Jaws
 - C. Insulation
 - D. Drive Rod/Threads
4. Verify jaws are closed; gently close if applicable. **Refer to picture 2**
5. Verify the handle is in the closed position and insert the threaded end of the single-use tip into the distal end of the shaft. **Refer to picture 3**
Note: Failure to completely close the jaws and/or handle prior to assembly may compromise the cutting/grasping performance.
6. While keeping handle closed, rotate the tip clockwise until it is secure. No gap should remain between the tip and shaft. **Refer to picture 4**
Note: If applicable, verify the shaft O-ring is covered and tighten further if not covered. No gap should remain between the tip and shaft. Failure to properly assemble the device may result in instrument malfunction and possible patient or user injury.
7. Remove tip protector from the end of the tip. **Refer to picture 5**
8. Open the tip to full aperture. If the jaws are difficult or not possible to move, remove the tip by rotating it counterclockwise and repeat Steps 4 through 6 above. **Refer to picture 6**
9. Keep the contact surface of the instrument clean during the operation. Wipe off any dried residue.
10. The tip mechanism is designed to open and close smoothly. Care should be taken not to forcibly open the jaws wider as undue stress may damage the insulation, or jaw components.

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Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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