





GENICON SINGLE USE CLIP APPLIER INSTRUCTIONS FOR USE

Including sterile packaged Single-Use 10 mm Multi-Fire Clip Applier

 **BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

IMPORTANT!

1. This booklet is designed to assist using this product. It is not a reference to surgical techniques.
2. This device was designed, tested and manufactured for single use patient use only.
3. Reuse of this device may lead to its failure and subsequent patient injury.
4. Reprocessing and/or resterilization of this device may create risk of contamination and patient infection.
5. Do not reuse, reprocess or resterilize this device.

DESCRIPTION

The GENICON A Ligator is a sterile packaged single-use clip applier designed to provide a means of ligation through an appropriately sized cannula. The device includes a rotatable shaft containing 20 titanium clips.

INDICATIONS FOR USE –

The GENICON Single-Use Clip Applier is intended for use on tubular structures or vessel wherever a metal ligating clip is indicated.

CONTRAINDICATIONS

1. The GENICON Laparoscopic single use clip appliers are NOT intended for contraceptive tubal occlusion.
2. The GENICON Laparoscopic single use clip appliers are NOT intended to be used on tissue structures or vessels upon which metal ligating clip would not normally be used.
3. These devices are intended for use only as indicated.

WARNINGS AND PRECAUTIONS

1. The use of GENICON laparoscopic and endoscopic products is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason. Contraindications relevant to individual GENICON products are noted in the specific sections.
2. Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions, and corrective actions in the event of failure.
3. This device is provided STERILE and is intended for SINGLE procedure ONLY. DISCARD AFTER USE. DO NOT RESTERILIZE.
4. Do not use any pre-sterilized products if the packaging sterile barrier is damaged.
5. Single Use devices may not be repaired, modified or re-processed.
6. Handle with care and store all products safely at temperatures less than 54 ° C and relative humidity between 0% and 80%.
7. Ensure that the clip is the correct size for the vessel or tubular structure being ligated.
8. Do not insert the Clip Applier through a trocar if a clip is present in the jaws. This may result in clip malformation, dislodged clips, or damage to the instrument. If a clip is present, fully squeeze the trigger against the handle, then fully release the trigger to release the clip from the jaws.
9. Prior to loading a clip in the jaws, ensure that the jaws are past the end of the trocar cannula.
10. Always check to see if a clip is seated in the jaw prior to firing!
11. Inspect the jaw tips to ensure the clip is fully advanced in the jaws.
12. Ensure that each clip is securely and completely positioned around the tissue being ligated.
13. Do not excessively twist or torque the instrument jaws when positioning or firing the instrument on a tubular structure or vessel. Excessive twisting or torquing may result in clip malformation.
14. Do not excessively apply a side load to the jaws that would cause them to partially collapse and potentially result in a clip malformation.
15. Excessive tissue manipulation with clip in jaws may result in clip dislodgement.
16. The jaws should be fully open and parallel upon initiating the firing of the instrument.
17. Do not push with excessive force on the proximal end of the instrument. Excessive force may result in clip malformation.
18. The trigger must be fully squeezed against the handle to ensure complete clip formation.
19. Ensure full release of the trigger after firing. A partial release of the trigger may disrupt clip feeding sequence and may result in clip malformation.
20. Avoid firing the instrument over another clip or instrument. Firing the instrument in this manner may distort or yield the instrument jaws.
21. If a clip is dislodged prematurely from the jaw tips or a clip fails to advance, remove the jaws from the targeted structure and fully squeeze and release the trigger to reset the device.
22. If a clip is present in the jaws and the clip applier needs to be removed from the patient, fully squeeze and hold the trigger against the handle while removing the clip applier through the trocar with the jaws in the closed position. Once the clip applier is removed, fully release the trigger to release the clip from the jaws. The Clip Applier is then ready for the next clip application.

POTENTIAL COMPLICATIONS FOR ALL GENICON® PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage/perforation, hemorrhage, hematoma, trocar site hernia, and septicemia/infection.

INSTRUCTIONS FOR USE

1. Open the package by using standard sterile technique.
2. Inspect the unit for defects – do not use if unit is defective. **Refer to pictures 1**
 - A. Jaws
 - B. Shaft
 - C. Rotation Knob
 - D. Trigger
 - E. Handle
3. Insert the Clip Applier (with NO clip loaded in the jaws) through an appropriately-sized cannula. The empty jaws will passively collapse as they pass through a cannula and reopen when completely through the cannula.
4. To load the first clip in the jaws: fully squeeze and release the trigger to reset the device. **Refer to pictures 2**
5. Load a clip into the jaws by slightly squeezing the trigger. **Refer to pictures 3**

NOTE: Inspect the jaw tips to ensure the clip is fully advanced in the jaws.
6. Position the jaws with the clip completely around the tubular structure or vessel to be ligated. **Refer to pictures 4**
7. The shaft can be rotated 360 degrees to facilitate visualization and accurate placement. **Refer to pictures 5**
8. Fully squeeze the trigger until it stops against the handle to completely form the clip on the targeted structure or vessel. **Refer to pictures 6**
9. After firing, fully release the trigger. A clip will not be loaded in the jaws until the trigger is squeezed again.
10. Check to ensure that each clip has been securely placed around the tissue being ligated.
11. The clip applier can be used to secure a catheter for cholangiography. During closure on the cystic duct and catheter, release the trigger prior to the trigger stopping against the handle. If the clip is dislodged from the jaw remove the jaws from the targeted structure and fully squeeze and release the trigger to reset the device.
12. To remove the instrument, ensure there is no clip remaining in the jaws, and withdraw the instrument from the cannula.

NOTE: The instrument contains a last clip lockout feature designed to increase the force required to close the trigger.

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This product may be covered by one or more U.S. patents and their foreign counter parts. Other Patents may be pending.
GENICON and the GENICON logo design are U.S. registered trademarks of GENICON Rights Reserved.

Disclaimer of Warranties on Re-Processing or Re-Sterilizing of Vendor's Disposable (Single-Patient) Products: Any reprocessing or re-sterilizing of disposable single patient devices voids Vendor's warranty. Vendor assumes no liability for reprocessed or re-sterilized products. Vendor warrants that reasonable care has been used in the manufacture of these disposable single patient devices and that they are free from defects in workmanship or materials at the time of shipment from Vendor. Vendor's sole obligation shall therefore be to repair or replace any device which it determines was defective at the time of shipment. Because no product is completely effective under all circumstances, and because the actual use and handling of these devices are beyond Vendor's control, Vendor cannot warrant for a good effect or against a bad effect in the application and use of these disposable single patient devices. The Purchaser therefore assumes all liability resulting from re-sterilization or re-processing of these disposable single patient products. Vendor therefore gives no warranty of merchantability or fitness for a particular purpose. Vendor will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to a disposable product that has been re-sterilized or re-processed. Vendor shall not be liable for incidental or consequential loss, damage or expense resulting from the use or application of these said products. This warranty is in lieu of all other warranties, whether implied, express, oral or written, and no individual or entity has the authority to vary the terms of this warranty.

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