



New Ideas for Modern Health Care

Comfort, Control and an Easy Mode of Entry



Hasson Port System

Access Ports: Hasson Port System

The GENICON Access Port Systems are designed to offer unrivaled versatility, efficiency, ease of use, and low cost.



Hasson with height-adjustable, angle suture wing with single-hand adjustable sliding function.

Industry first ergonomic palm grip for a more comfortable insertion



Soft and malleable silicone cone conforms to skin and abdominal tissue



Multiple versatile suture slots

Hasson / Blunt Port System

| | Diameter | Length | Qty/Box | Part Number |
|---------------------|----------|--------|---------|-------------|
| With Smooth Cannula | 12mm | 100mm | 5 | 400-012-131 |